


Framework Management Digital



AMANNGIRRBACH

 **ceramill map**

 **ceramill mind**

 **ceramill motion**

 **ceramill m-plant**

 **ceramill artex®**

 **ceramill m-center**



AMANNGIRRBACH

MODEL MANAGEMENT



Registration | Articulation

Model fabrication

Dosing | Mixing | Cleaning

FRAMEWORK MANAGEMENT



DIGITAL

Outsource Digital

Inhouse Digital

ANALOG

Inhouse Manual



FRAMEWORK MANAGEMENT

 **ceramill mall**



**Outsource
Digital**




**Inhouse
Digital**



**Inhouse
Manual**



Looking for simplicity,
precision and affordability?
We understand your needs!



In-house CAD-CAM at
affordable prices



Simple and user-friendly
operation



Precision from start to finish for the perfect fit

FRAMEWORK MANAGEMENT DIGITAL

OVERVIEW		F 6	
SYSTEM COMPONENTS	Ceramill Map100	F 12	
	Ceramill Map300	F 13	
	Ceramill Mind	F 14	
	Ceramill M-Plant	F 17	
	Ceramill TI-Connect	F 19	
	Ceramill Artex	F 15	
	Ceramill Upload-Tool	F 21	
	Ceramill ZI	F 27	
ORDER INFORMATIONS		B 29	

KNOW-HOW

	Ceramill Helpdesk	KH 36	
	Training courses	KH 37	
	Webinar	KH 39	
	D-Lab24	KH 40	

SERVICE

	Customer service	S 42	
	Technical service	S 42	
	Service partner	S 42	
	AG Live Labs	S 43	
	General Information	A 43	

OVERVIEW

ceramill mall

Ceramill Mall
system components:

ceramill map100

The model scanner with a perfect price performance ratio



ceramill map300

The high performance model scanner with connection to articulator



ceramill artex®


The interface between manual and digital dental technology



<p>_Economical introductory model</p>	<p>_Artex® CR synchable with virtual articulator</p>	<p>_Range of functions 1:1 to the actual Artex® CR</p>
<p>_Functional operation concept</p>	<p>_Fully automatic scanning of bridges up to 14 units</p>	<p>_Perfect synchronisation between the actual and virtual articulator</p>
<p>_Bite registration, gingiva and wax-up scans possible with ease</p>	<p>_Unlimited indication width</p>	<p>_Software-based 1:1 visualisation of the Artex® CR (incl. excursion simulation)</p>




OVERVIEW



ceramill mind

The brain of the system:
The construction software




ceramill motion

Manageable, versatile and accurate:
The CAM milling machine



ceramill m-center

Reliable and precise:
The milling center



_ Economizing time through intuitive dental technical workflow

_ Large functional scope

_ Automatic preparation margin recognition
_ Comprehensive indication range

_ Amortisation in record time

_ Wide spectrum of indications regarding materials and size

_ Superior CNC technology for high productivity

_ Completion of indication range without additional investment and stock keeping

_ Quicker through-put times with high quality & precision

_ Safer and simpler flow of work

OVERVIEW

outsources digital


In-house construction,
external milling




METHOD



Ceramill Cockpit:
digital control panel for the overall Ceramill production process, selection of indication mode.




Scanning of the model via high precision strip-light projection.



- 1 Automatic recognition of the preparation margin.
- 2 Calculation of the anatomically reduced framework based on the full anatomic construction.

At the click of a button the saved construction data are sent to the Ceramill M-Center.



After the data input control, the construction is assigned, according to Material, to the designated manufacturing technique and the production process is calculated.



The construction is then produced either by milling with a CNC milling machine or using the laser melting technique.



After the quality control the framework is packaged and dispatched.

SYSTEMS AND PRODUCTS



ceramill map100

F 12


- _ Low-cost entry due to partially automatic features as opposed to fully automatic
- _ Easy re-scanning possible
- _ Main application: Small pieces of work (up to 6 units, larger pieces of work are also possible)



ceramill map300

F 13


- _ Fully automatic scanner, also registers bridges up to 14 units quickly and easily
- _ Compatible with the virtual articulator Ceramill Artex® CR for automatic, full anatomical framework design.



ceramill mind

F 14

- _ Accurate and quick recognition of the preparation margin
- _ Intuitive dental technical work-flow which provides easy operation



ceramill zi units

F 25

- _ Zirconium oxide milled and sintered in the Ceramill M-Center



ceramill np units

F 25

- _ CoCr bridges produced using laser-melting techniques in the Ceramill M-Center

d-lab24.com

F KH 40

Digital instructions, videos, forums and much more training material on AG methods.



OVERVIEW

inhouse digital

In-house digital construction and milling




METHOD



Ceramill Cockpit: digital control panel for the overall Ceramill production process, selection of indication mode.



Scanning of the model via high precision strip-light projection.



- ① Automatic recognition of the preparation margin.
- ② Calculation of the anatomically reduced framework based on the full anatomic construction.



Placing the construction in the desired blank, adaptation of bars and calculation of milling paths.



Milling the construction in the milling machine Ceramill Motion



Stained and sintered bridge

SYSTEMS AND PRODUCTS



ceramill map100

F 12

- _ Low-cost entry due to partially automatic features as opposed to fully automatic
- _ Easy re-scanning possible
- _ Main application: Small pieces of work (up to 6 units, larger pieces of work are also possible)



ceramill map300

F 13


- _ Fully automatic scanner, also registers bridges up to 14 units quickly and easily
- _ Compatible with the virtual articulator Artex® CR for automatic, full anatomical framework design.



ceramill mind

F 14

- _ Accurate and quick recognition of the preparation margin
- _ Intuitive dental technical work-flow which provides easy operation



ceramill match

F 24

- _ Quick and easy positioning of the units in the blank - efficient use of the blank



ceramill motion

F 22

- _ For dry finishing zirconium oxide and acrylic
- _ Mills up to 14 unit bridges
- _ Optimised milling process, lower milling times, high productivity

d-lab24.com

KH 40

Digital instructions, videos, forums and much more training material on AG methods.



OVERVIEW

ceramill material

Overview of the wide range of materials and applications covered by Ceramill Mall:

Ceramill Framework units:	The most important features at a glance:	The most important indications:	
Ceramill ZI units (zirconium oxide) 	<ul style="list-style-type: none"> _ High strength, rigidity, and biocompatibility _ On request tooth coloured framework by staining _ Conventional and adhesive fixation 	<ul style="list-style-type: none"> _ Crowns and bridges in the anterior and posterior region _ Telescope and conical crowns _ Abutments and attachments 	
Ceramill COMP units (Composite) 	<ul style="list-style-type: none"> _ Nano composite _ High flexural strength _ Metal-free _ Ready-made material from industrial production 	<ul style="list-style-type: none"> _ Permanent prostheses: up to 3 unit bridges _ Suitable for veneer bonding <p>Approved for up to: 14-unit bridges, 3 pontics</p>	
Ceramill TEMP units (dyed PMMA) 	<ul style="list-style-type: none"> _ Dyed resin _ Industrially pre-fabricated material 	<ul style="list-style-type: none"> _ Temporary appliances _ 2 intermediate elements _ Suitable for veneer bonding 	
Ceramill PMMA units (acrylic which burns out without residue) 	<ul style="list-style-type: none"> _ Transparent acrylic _ Burns out without leaving a residue _ Industrially pre-fabricated material 	<ul style="list-style-type: none"> _ Framework for intraoral try-in _ Framework for casting technique and press-on technique 	
Ceramill WAX units 	<ul style="list-style-type: none"> _ Burns out without leaving a residue _ Industrially pre-fabricated material 	<ul style="list-style-type: none"> _ Framework for casting technique and press-on technique 	



OVERVIEW



Ceramill Framework units:	The most important features at a glance:	The most important indications:	
<p>Ceramill GCER LS 2 units (glass ceramic)</p> 	<ul style="list-style-type: none"> _ High aesthetics _ High strength _ Industrially pre-fabricated material _ For conventional and adhesive fixation 	<ul style="list-style-type: none"> _ Inlays, Onlays, Veneers _ Single crowns 	
<p>Ceramill NP L units (CoCr - Lasergesintert)</p> 	<ul style="list-style-type: none"> _ High strength, rigidity _ Complex frameworks possible _ Beryllium and nickel-free _ For veneering with conventional metal bonding porcelains 	<ul style="list-style-type: none"> _ Crowns and bridges in the anterior and posterior region _ Telescope and conical crowns _ Approved as a framework for veneering up to 7 units 	
<p>Ceramill NP M units (CoCr - Geprüft)</p> 	<ul style="list-style-type: none"> _ High strength, rigidity _ Complex frameworks possible _ Beryllium and nickel-free _ For veneering with conventional metal bonding porcelains 	<ul style="list-style-type: none"> _ Crowns and bridges in the anterior and posterior region _ Telescope and conical crowns _ Frameworks for veneering or fully anatomical 	
<p>Ceramill TI alloy units (TiAlNb)</p> 	<ul style="list-style-type: none"> _ High strength, rigidity, and biocompatibility _ Homogenous structure _ Industrially pre-fabricated material, no casting defects 	<ul style="list-style-type: none"> _ Crowns and bridges in the anterior and posterior region _ Telescope and conical crowns _ Individual abutments on titanium bases 	
<p>Ceramill TI units (Grade 2)</p> 	<ul style="list-style-type: none"> _ High strength, rigidity, and biocompatibility _ Homogenous structure _ Industrially pre-fabricated material, no casting defects 	<ul style="list-style-type: none"> _ Crowns and bridges in the anterior and posterior region _ Telescope and conical crowns _ Individual abutments on titanium bases 	

SYSTEM COMPONENTS



ceramill map100

The part-automatic entry-level scanner - uncompromising scanning quality

All scanners in the Ceramill Map range are characterised by their high resolution data, which is generated by the strip light projection. Highly sensitive 3D sensors ensure of highly accurate images of the model (< 20 µm).

Another standard quality is the functional range of the background software: Apart from scanning the model, it is also capable of scanning the bite registrations, gingiva and wax-ups.

The low-cost basic scanner Ceramill Map100 is ideally equipped for 90% of all prosthetic restorations (6 units). Larger restorations are also possible.

The scanning field comprises 50x36x40 mm. The model is positioned by hand. The end position can be checked live on the monitor.



- _The automatic axis has been replaced by a manually manoeuvrable axis, thus providing the low entry-level price
- _Individual scanning strategy (manual model alignment), saves time
- _The model can be adjusted by hand in the scanner. Re-scanning takes place with just one mouse click
- _The model, the live image on the monitor and the constructed scan all move synchronously so that the missing areas can be directly re-scanned
- _Scanning can take place while the door is open due to the strip light projection
- _Main application: 1-6 unit bridges, 14 unit bridges also possible
- _Bite registration, gingiva and wax-up scan possible for optimised framework fabrication
- _The scanner comes with an open interface, scans (stl-files) can also be downloaded into other CAD software



The model can be adjusted in the scanner by hand. Side effect: Waiting times during the scanning process are reduced.



When the model is turned in the scanner, the live image and the scan also turn automatically on the monitor.



Scan and live image - easy re-scanning via mouse-click or button on the scanner.



SYSTEM COMPONENTS



ceramill map300

The fully automatic scanner - uncompromising quality and handling

Apart from the same qualities as the Map100, the Ceramill Map300 is particularly characterised by its ease of use.

The scan field is identified and assessed automatically, which means that large spanned bridge frameworks can be quickly and efficiently documented.

Articulated models can also be easily processed by this type - one of the prerequisites when the "virtual articulator" comes into action with the construction software.




- _The fully automatic strip light scanner scans up to 14 unit bridges quickly and easily
- _Articulated models can be scanned in relation to the articulator, then downloaded into the CAD software in the Ceramill Artex® CR for the automatic fabrication of fully anatomical frameworks
- _Bite registration, situation model, gingiva and wax-up-scan possible for optimal framework fabrication
- _Automatic user-guidance through the scan programme for easy and safe operation
- _The scanner comes with an open interface, scans (stl-files) can also be downloaded into other CAD software



Articulated models in the Ceramill Fixator and Artex® CR. The models were synchronised using the Splitex® key.



Ceramill Fixator with articulated model in the Ceramill Map300 (Symbol illustration of the Ceramill Fixator).



Ceramill Map300 with model in the Ceramill Fixator - for lossless transfer of the model situation.

SYSTEM COMPONENTS



The intelligent Construction Software - The "spirit" of framework fabrication

A software which is designed to follow the dental technical work-flow is the main prerequisite for the greatest possible acceptance of a digital medium. The CAD software Ceramill Mind was developed in close collaboration with dental technicians - and its further development is confirmed. Its contents leave nothing to be desired: Precise recognition of the prepared margin, automatic bridge and connector design, tooth library compilation and the open system structure are just some of the highlights.



Simultaneous designing in the upper and lower jaw

- _Large indication spectrum (fully anatomical and anatomically reduced crowns and bridges, inlays/onlays, reducing wax-ups, press-over, telescopes, virtual articulator)
- _Simultaneous designing in the upper and lower jaw
- _Library teeth automatically adapt to the scanned diagnostic model
- _Automatic defining of the preparation margins
- _Intuitive dental technical work-flow for comfortable and reliable use
- _The virtual articulator simulates the mandibular excursions and automatically constructs a full anatomical framework proposal, according to the **dynamic** occlusion, therefore reducing the necessity to grind the surfaces after milling
- _The order button integrated in the programme makes sending the construction data simple
- _Open for every type of construction data (File type: stl), increases flexibility




COURSES


Please also note our course offers on this subject on Page 37



14-unit bridges with 5 individual abutments



Automatic display of the connector dimensions



Automatic adaptation of the library teeth to the teeth of the scanned diagnostic model



SYSTEM COMPONENTS



Virtual Artex® CR as Upgrade for Ceramill Map300 and Ceramill Mind.
The functional interface between manual and digital prosthetic dentistry

With manual production of dental prosthesis working with the articulator is standard for dental laboratories. In order to achieve the same quality of the works virtually, it is only logical and consistent to enable this by means of a CAD-CAM system. The virtual articulator "Ceramill Artex" serves as a bridge between manual and digital techniques:

The model pair in the Artex articulator is transferred to the Map300-scanner while holding the same Artex mode by means of the Ceramill transferkit; it is subsequently scanned-in in the appropriate proportion.

The motion mode and free zones of the Artex®CR are thus brought into line on the same modes - digitally and manually.

Disturbing structures may already be removed thus reducing time-consuming grinding in the patient's mouth to a minimum.



- _ Fully visualised Artex® CR for a quick introduction to the digital world
- _ The virtual articulator offers the same functional scope as compared to the real Artex® CR (Adjustment modes of the horizontal inclination of condylar guidance (Bennett Angle; Retrusion; Immediate Side Shift)
- _ The transfer of the models by means of the Ceramill® Fixator ensures the precision at the functional interface between manual and digital techniques.
- _ The calculation of the fully anatomical construction is dynamic and static under consideration of the antagonists and the adjusted values of the articulator.
- _ Space for the porcelain built-up is automatically foreseen during the construction - thus an optimal framework basis is established for veneers with high stability and a consistent layer thickness.




VIDEO

Demonstration video at www.amangirrbach.com



Models in the Ceramill-Fixator



Models in the Ceramill fixator as a transfer interface in the Ceramill Map300



Adjustment modes at the virtual Artex® CR

SYSTEM COMPONENTS



Step-by-Step

The Ceramill Artex® virtual articulator offers exactly the same setting options as the manual version (Artex®CR articulator). The articulator can be set in the same way as the original using a software mask. Adjustments to the articulator setting are completed onscreen and animated in real time on the Ceramill Artex®. This enables an immediate visual control of the settings and therefore makes the virtual articulator “functional”.




A) Condyle actual Artex®CR




B) Condyle virtual (0 degrees)



C) Condyle virtual (30 degrees)




A) Artex®CR condyle in centric position viewed from below




B) Virtual Artex®CR condyle in centric position as starting point for each excursion




C) Artex®CR condyle in (animated in real time) lateral excursion



A) The CAD design with visual marking of contact and penetration areas to the opposing model before use of the virtual articulator (calculation of the dynamics)



B) Ceramill Artex® in function - static reduction of the CAD design in the functional surfaces



C) The result of the dynamically automated operation of the Ceramill Artex®




SYSTEM COMPONENTS



ceramill m-plant

Abutment Module Upgrade for Ceramill Mind

Ceramill M-Plant is an upgrade module for the construction software Ceramill Mind and functions as a complement to construct individual hybrid abutments made from titanium and zirconium oxide. With the adhesion of the titanium basis and the individual abutment the placement of a permanent and stable abutment in the implant is guaranteed.



Data software is available for the following implant systems:

- ____ Ceramill ti-connect for Nobel Biocare Replace Select®
- ____ Ceramill ti-connect for Nobel Biocare Nobel Active™
- ____ Ceramill ti-connect for Nobel Biocare Brånemark®
- ____ Ceramill ti-connect for Biomet 3i Osseotite® Certain®
- ____ Ceramill ti-connect for Biomet 3i Osseotite®
- ____ Ceramill ti-connect for Straumann® BoneLevel
- ____ Ceramill ti-connect for Straumann SynOcta®
- ____ Ceramill ti-connect for Zimmer Tapered Screw-Vent®
- ____ Ceramill ti-connect for Astra Tech OsseoSpeed®
- ____ Ceramill ti-connect for Dentsply-Friadent Frialit/Xive®
- ____ Ceramill ti-connect for Camlog
- ____ Ceramill ti-connect for Prowital
- ____ Ceramill ti-connect for Thommen Medical

- ____ Construction of all components required for individualised implant prosthetics are possible in one step with this software: Abutment implant, anatomical framework, if necessary wax frameworks for the over press technique.
- ____ Using the hybrid technique (adhesive technique) a safe and durable fit of the abutment on the implant
- ____ Finishing of each individualised abutment using the Ceramill Motion or Ceramill M-Center
- ____ All produced from one source - software and titanium bases are perfectly matched and compatible with one another for save processing
- ____ Construction can be started immediately, even when the Ti bases have not yet been received in the laboratory
- ____ Emergence profile can be custom fitted to the gingiva
- ____ Ti-Connect Ti bases can also be used for Ceramill Base and Ceramill Multi-x
- ____ Biocompatible restoration, no black metal margins in the patient's mouth
- ____ Aesthetic advantages, the gingival emergence profile of the abutment can be customised for each case, perfect gingiva, perfect papillae after healing



SYSTEM COMPONENTS



4 steps to a finished customised abutment

1 Model preparation




Model with removable gingival mask and integrated laboratory implant.



Fit the scan body onto the laboratory implant.
Note: Check that the scan body fits exactly and fix it in position using a screw. Recommended: order a separate screw extra for each scan body for fixing in position and subsequent adhesive retention.

2 Scanning




The model is scanned in the Ceramill Map100 or Map300 scanner. Simply follow the instructions of the software.

3 Designing




Designing the abutment using the Ceramill M-Plant software. The individual design steps are illustrated in the video tutorial on the homepage www.ceramill-m-center.com. The abutment can be fabricated in the Ceramill M-Center or using the Ceramill Motion.

4 Adhesive retention



Protect the connection geometry of the titanium base (A) with separating agent and fix in position on the laboratory implant using a screw. Sandblast the adhesive surface of the base (B) and abutment using Al₂O₃ blasting medium 50 µm at max. 2 bar. Clean the adhesive surfaces thoroughly. Seal the screw channel with wax. RelyX™ Unicem® (3M Espe), Panavia® F2.0 (Kuraray) or other similar luting materials are recommended for adhesive retention (use a metal primer, if necessary). Adhere to the manufacturer's instructions.



Apply adhesive to the upper edge (B) of the titanium base.
Rotate the abutment when fitting it to the titanium base, to ensure distribution of the adhesive. Once resistance is felt, detect the final position by rotating the abutment.
The abutment must fit flush with the titanium base. Remove large amounts of adhesive residue immediately.



Remove the excess after the adhesive has cured using a silicone polisher. Remove the wax in the screw channel and clean the abutment together with the titanium base.



COURSE

Webinar available: Information can be found on Page 37

INFO

Comprehensive video tutorial with verbal instructions on www.m-center.com (only for registered customers)



SYSTEM COMPONENTS



ceramill ti-connect

Titanium base for all conventional implant systems

Range	Range one Ceramill ti-connect for Astra Tech® OsseoSpeed®				Range two Ceramill ti-connect for Biomet 3i® Osseotite® Certain®			
	Kit a 3,5/4,0 mm	Kit b 4,5/5,0 mm			Kit a 3,4 mm	Kit b 4,1 mm	Kit c 5,0 mm	
Titanium base incl. screw								
	792111	792112			792211	792212	792213	
Scan body								
	792121	792122			792221	792222	792222	
Laboratory implant								
	792131	792132			792231	792232	792233	
Screw								
	792141	792142			792241	792241	792241	
Range	Range three Ceramill ti-connect for Straumann® BoneLevel®				Range four Ceramill ti-connect for Nobel Biocare® Nobel Active™			
	Kit a 3,3 mm	Kit b 4,1/4,8 mm			Kit a 3,5 mm	Kit b 4,3/5,0 mm		
Titanium base incl. screw								
	792311	792312			792411	792412		
Scan body								
	792321	792322			792421	792422		
Laboratory implant								
	792331	792332			792431	792432		
Screw								
	792341	792341			792441	792442		
Range	Range five Ceramill ti-connect for Straumann® SynOcta®				Range six Ceramill ti-connect for Nobel Biocare® Replace Select®			
	Kit a 3,5 mm	Kit b 4,8 mm	Kit c 6,5 mm		Kit a 3,5 mm	Kit b 4,3 mm	Kit c 5,0 mm	Kit d 6,0 mm
Titanium base incl. screw								
	792511	792512	792513		792611	792612	792613	792614
Scan body								
	792521	792522	792523		792621	792622	792623	792624
Laboratory implant								
	792531	792532	792533		792631	792632	792633	792634
Screw								
	792541	792542	792542		792641	792642	792642	792642

Frailiti/Alive® is a Registered Trademark of Dentsply Friadent GmbH or affiliated companies. | OsseoSpeed® is a Registered Trademark of Astra Tech AB or affiliated companies. | Tapered Screw-Vene® is a Registered Trademark of Zimmer Dental GmbH or affiliated companies. | Osseotite® and Osseotite® Certain® are Registered Trademarks of Biomet 3i, LLC or affiliated companies. | Nobel Active® and Replace Select® are Registered Trademarks of Nobel Biocare Management AG or affiliated companies. | BoneLevel® and SynOcta® are Registered Trademarks of Straumann AG or affiliated companies.

All original implants can be obtained from the respective manufacturers. The titanium bases/interfaces should also be ordered directly from the respective manufacturers for Thommen®, Camlog® and Prowital®.



SYSTEM COMPONENTS



Range	Range seven Ceramill ti-connect for Biomet 3i Osseotite®			Range eight Ceramill ti-connect for Nobel Biocare Brånemark®			
	Kit a 3,4 mm	Kit b 4,1 mm	Kit c 5,0 mm	Kit a 3,5 mm	Kit b 4,1 mm	Kit c 5,1 mm	
Titanium base incl. screw	 792711	 792712	 792713	 792811	 792812	 792813	
Scan body	 792721	 792722	 792722	 792821	 792822	 792823	
Laboratory implant	 792731	 792732	 792733	 792831	 792832	 792833	
Screw	 792741	 792741	 792741	 792841	 792842	 792843	
Range	Range nine Ceramill ti-connect for Zimmer Tapered Screw-Vent®			Range ten Ceramill ti-connect for Dentsply Friadent Frialit/Xive®			
	Kit a 3,5 mm	Kit b 4,5 mm	Kit c 5,7 mm	Kit a 3,4 mm	Kit b 3,8 mm	Kit c 4,5 mm	Kit d 5,5 mm
Titanium base incl. screw	 792911	 792912	 792913	 7921011	 7921012	 7921013	 7921014
Scan body	 792921	 792922	 792923	 7921021	 7921022	 7921023	 7921023
Laboratory implant	 792931	 792932	 792933	 7921031	 7921032	 7921033	 7921034
Screw	 792941	 792941	 792941	 7921041	 7921041	 7921041	 7921041

Frialit/Xive® is a Registered Trademark of Dentsply Friadent GmbH or affiliated companies. | Osseo Speed® is a Registered Trademark of Astra Tech AB or affiliated companies. | Tapered Screw-Vent® is a Registered Trademark of Zimmer Dental GmbH or affiliated companies. | Osseotite® and Osseotite® Certain® are Registered Trademarks of Biomet 3i, LLC or affiliated companies. | Nobel Active® and Replace Select® are Registered Trademarks of Nobel Biocare Management AG or affiliated companies. | BoneLevel® and SynOcta® are Registered Trademarks of Straumann AG or affiliated companies.




SYSTEM COMPONENTS



ceramill upload-tool

The convenient interface to the Ceramill Motion and Ceramill M-Center for 3-Shape and Dental Wings users

The "Ceramill Upload Tool" software module allows users of open scanners (3Shape®, Dental Wings®) to fabricated design data in their own laboratory using the Ceramill Motion milling machine or in the Ceramill M-Center manufacturing centre.



A requirement for this is the compatibility of external data.

Simply test in three easy stages whether your 3Shape or Dental Wings' scanner data is compatible with the Ceramill Motion and Ceramill M-Center.



Simply complete the online form (at www.m-center.com) and send us your open STL file (ZIP compressed*).




Your data will be subjected to a compatibility check for the Ceramill Motion as well as the Ceramill M-Center.



After the successful test our field service will contact you (**within 2 working days** following the day of your data dispatch).



Simply complete the form for the test data upload of STL data.



Load your STL data at the press of a button into the Ceramill Upload-Tool.



Milling of the digital data in the Ceramill Motion or Ceramill M-Center.

* ZIP-packed = STL-files packed with WinZip®-packing program - Available at no charge on the Internet

SYSTEM COMPONENTS



ceramill motion

The compact milling unit with intelligent 4-axis control - for rapid and accurately fitting framework fabrication in zirconia, resin, non-precious metal and wax


The compact CNC* milling machine Ceramill Motion gets things moving.

Clean, efficient and precise - this CAM unit mills most framework situations in the laboratory.

Using Ceramill Match, the works aligned with the blank are converted into milling data at lightning speed.

The compact external dimensions and an appealing design make the Ceramill Motion an attractive looking work horse to have in the laboratory.

The 4-axes compact milling unit guarantees a highly precise fit, even with undercut areas.



- _ For dry milling of zirconia, resin, non-precious metal and wax
- _ Mills up to 14 unit bridges
- _ Optimised milling process, therefore lower milling times, high productivity
- _ Small amount of space required, fits into every laboratory
- _ Tool holding fixture with automatic tool changer
- _ Open interface for an open scanner (3Shape®, Dental Wings®)
- _ Quick and easy positioning of the units in the blank - efficient use of the blank
- _ Highly precise Jäger spindle with a true-running accuracy of < 0,004 mm
- _ Smooth working possible by means of integrated monitoring electronics of the axis position
- _ Air jet function and tool cooling by integrated air nozzles at the spindle for a long operating time
- _ Automatic tool length measurement and broken tool detection
- _ Reserve tools: If the first tool set is worn, the machine automatically changes to a new set
- _ Minimum maintenance required



SYSTEM COMPONENTS



ceramill motion

Presence sensor
for milling tools

Tool-length sensor
incl. fracture check and calibration

Tool magazine
with automatic tool changer

Blank holder

Interior lighting
for visual control of milling

Removable extraction tray
for easy cleaning


Connection for extraction protects the internal mechanism against contamination
Optimal for use with the Ceramill Airstream, but can also be used for central extraction





Function assignment:
- Reference run
- Start/Pause programme
- Workpiece control
- Extraction

Multifunction button


Function assignment:
- Tool holder
- Service position
- Spindle ventilation



Engage the 4th axis to any position



Milling of undercuts possible



SYSTEM COMPONENTS




ceramill motion

Inhouse milling with premium performance, usability and precision



Only the best is good enough. The centrepiece of each Ceramill Motion is a high frequency spindle of the worldwide known German company Alfred Jäger that - with 40 years of experience in the area of engine spindles and application technology - ensures profound expert knowledge and reliability to meet the extraordinary high performance requirement of Ceramill Motion.



ceramill match

The automatic operator guidance and the transparent user interface of the Ceramill Match CAM software form the basis for a reliable and easy operation. No CAM or milling know-how is required to use it. Even users with little experience may quickly and easily establish the milling programs to manufacture crowns and bridge frameworks. An elaborated collision control (and evasion) of Ceramill Match® ensures a high degree of process reliability.



Tool magazine with automatic tool changer incl. reserve tool places.



High precision in zirconia, resin, non-precious metal and wax.



SYSTEM COMPONENTS



ceramill m-center


The manufacturing centre for all indications and materials:
Versatile, reliable and fast

The prosthetic restoration still bears the individual technician's signature - the milling center will only produce a piece of work which has been previously approved by him.

A procedure using state-of-the-art technology and specially trained processing engineers. Not without dental technical quality control, which takes place before the finished framework is dispatched.

Naturally, various AG support media is available on all aspects in the Ceramill M-Center: Online help, tracking the status of each order, the Ceramill help desk.

www.ceramill-m-center.com



- _The wide spectrum of different materials available (zirconium oxide, glass ceramic, CrCo, titanium, composite, PMMA) cover almost the entire framework fabrication spectrum
- _Diverse indications (anatomically reduced and fully anatomical crowns and bridges, telescope crowns, inlays/onlays, press-on technique, individual abutments)
- _Quick processing times - good planning safety
- _State-of-the-art production techniques produce optimal and precise frameworks
- _Video-tutorials, CAD-CAM forum and exclusive services for M-Center customers



Ceramill M-Center is equipped with the latest lasermelting and ultrasonic[®] milling technology and produces to the highest quality standards.

Quality control by professionals.

SYSTEM COMPONENTS



ceramill-m-center.com

Your interface to the AmannGirrbach manufacturing centre

The Ceramill M-Center Website is designed in particular for customers with an AmannGirrbach CAD/CAM system. It provides the customer with basic information on the product, the method and the manufacturing process. For registered customers, it offers a wide range of additional services like transfer of know-how and order processing. Hardware and software registration, user support, monitoring the order status (*no ordering function) as well as package tracking guarantee customers full control of the digital outsourcing process chain.



One-click upload of construction- and order data using the upload-button of the Ceramill mind software. Therefore there is no need of any additional order function via the website.

Only for registered customers in the download area.




Order processing



Download area (only for registered customers)

In the following languages:

- German
- English
- French
- Italian
- Spanish



CAD-CAM user forum (only for registered customers)



Information material/indications



SYSTEM COMPONENTS



ceramill zi

Presintered Y-TZP zirconium-oxide blanks for machining without risk of splintering or damaged edges


Ceramill zirconia blanks already undergo strict test procedures according to fixed parameters during the manufacturing process in our in-house production facility. In this way AmannGirrbach ensures the permanent high-quality standard of its in-house production.

The pre-sintered blanks can be perfectly machined - they do not splinter and offer excellent edge stability. All blanks are marked in batches with their individual enlargement factors, which are later transferred to the milling unit Ceramill.




- ____ TOSOH Powder
- ____ All Biaxial Pressed
- ____ 1,300 MPa Strength
- ____ FDA & Health Canada Approved
- ____ IdentCERAM Stickers


- ____ High-grade certified raw material
- ____ Extremely high bending strength, optimal edge stability and machinability
- ____ Due to their reusability, the blanks can be used efficiently so that unmachined areas can be processed later
- ____ Optimal fitting due to encoded indication of the enlargement factor on the blanks
- ____ Individual coloring due to 4 different degrees of brightness of the Ceramill Liquid dyeing solutions
- ____ Can be used with all usual zircon veneering ceramics
- ____ Designed exactly according to the course of the dental arch
- ____ The ideal blank for large spanned bridges
- ____ Covers up to 90% of all work in a laboratory, with no material wastage
- ____ Different heights




Possible unit per Ceramill Motion blank



At approx. 1300 Mpa flexural strength, Ceramill ZI lies in direct competition with the leading specialists.



Ceramill ZI has an excellent Weibull-modulus of 14.2 and therefore proves a high degree of reliability during use.



Flexural strength of Ceramill ZI non-coloured and coloured with Ceramill liquid; AmannGirrbach, 2009

SYSTEM COMPONENTS



A fully automatic, high performance furnace for the final sintering of stress-free frameworks

The Ceramill Therm is a high temperature furnace with additional features.

The milled Ceramill zirconium-oxide frameworks are dense-sintered with the Ceramill Therm and thus obtain their final density and the resulting excellent material properties. For sintering, the objects are placed onto sintering beads, which ensures a frictionless sintering process and thus distortion-free frameworks. The Ceramill Therm offers maximum process reliability due to constant temperature control, even temperature distribution in the firing chamber and notification in case of termination of the sintering programme due to e.g. power failure. As a result of this, the user is able to safely control if the final density and thus strength of the frameworks has been achieved.

The user has 4 sintering programme locations at his disposal; one of them can be individually programmed.



- _Maximum process reliability due to constant temperature control, even temperature distribution in the firing chamber and notification in case of termination of the sintering programme
- _Maximum process reliability due to optimally coordinated, fully-automated sintering programmes for different restoration sizes
- _4 sintering programme locations; one of them individually programmable by the user
- _2 stackable sintering bowls for maximum utilisation of the furnace
- _Minimum required space and installation time (supply required)



For maximum utilisation of the furnace, the Ceramill Therm exhibits a large firing chamber, in which 2 sintering bowls can be stacked on each other and thus the double amount of frameworks can be sintered in one cycle. With the aid of the sinter forceps, the sinter shell can be transferred easily and safely in and out of the furnace.





ORDER INFORMATIONS



Ceramill Map100

Ceramill Mind

Ceramill Outsource Digital 100 | SET

- 179003 - Ceramill Map100 (part automatic scanner)
 - Ceramill Mind (CAD Software)
 - PC incl. monitor, keyboard, mouse, virus scan software
- without furnace and extraction



Ceramill Map300

Ceramill Mind

Ceramill Outsource Digital 300 | SET

- 179004 - Ceramill Map300 (full automatic scanner) incl. Ceramill Transferkit
 - Ceramill Mind (CAD Software)
 - Ceramill Artex (virtual Articulator, upgrade for Ceramill Mind)
 - PC incl. monitor, keyboard, mouse, virus scan software
- without furnace and extraction



Ceramill Mind

Ceramill Motion

Ceramill Inhouse Digital Motion | SET

- 179005 - Ceramill Motion incl. Ceramill Match (Milling device incl. CAM Software)
 - PC incl. monitor, keyboard, mouse, virus scan software
- without furnace and extraction



Ceramill Map100

Ceramill Mind

Ceramill Motion

Ceramill Inhouse Digital 100 | SET

- 179001 - Ceramill Map100 (part automatic scanner)
 - Ceramill Mind (CAD Software)
 - Ceramill Motion incl. Ceramill Match (Milling device incl. CAM SW)
 - PC incl. monitor, keyboard, mouse, virus scan software
- without furnace and extraction



Ceramill Map300

Ceramill Mind

Ceramill Motion

Ceramill Inhouse Digital 300 | SET

- 179002 - Ceramill Map300 (full automatic scanner)
 - Ceramill Mind (CAD SW)
 - Ceramill Artex (virtual Articulator, upgrade for Ceramill Mind)
 - Ceramill Motion incl. Ceramill Match (Milling device incl. CAM Software)
 - PC incl. monitor, keyboard, mouse, virus scan software incl. Ceramill Transferkit
- without furnace and extraction





Technical Data:



Ceramill	Map 100	Map300	Motion
Art.Nr.	179100	179110	179200
Dimensions D/W/H (mm)	426 x 328 x 603	426 x 414 x 720	465 x 430 x 730
Weight (kg)	33	57	62
Power supply (V/A)	230/1,25	230/1,25	230/3,15 100; 115/6,3
E-fuse	T1, 25A	T1, 25A	
Output (W)	50	50	250
Compressed air connection			dry, clean compressed air, 6 bar max. 50 L/min
Engine speed (U/min ⁻¹)			60.000
Torque (Ncm)			4
Collet chuck (Ø mm)			3
Sound level (max. dbA)			60
Accuracy (µm)	<20	<20	<10
Axes	2	4	4
Recommended installation site	Table, no direct sunlight	Table, no direct sunlight	Table, no direct sunlight
Recommended temperature	18-30 °C	18-30 °C	18-30 °C

System requirements:


- _ Internet connection (minimum DSL) to guarantee data transmission to the M-Center and remote maintenance of the system
- _ Flat rate is recommended
- _ A network cable for connecting the router / modem and PC.
- _ An Internet connection must be available on the day of installation. If the customer does not know how to create an internet connection, it must be ensured that an Internet specialist is on site on the day of installation. AG will not connect the system to an already existing network at the customer.
- _ The customer must guarantee maintenance and any problem solving relating to the network / Internet connection.



Technical Data PC


Processor speed/model:	Intel Xeon Quad 2.67 GHz
Memory:	4 GB RAM
Hard drive:	232 GB
Graphics:	T1, 25A
Operating system:	Windows XP, Recovery Software, Virus protection MC Afee
Display:	19 Zoll
Resolution pixel:	1280x1024
Recommended installation site:	Away from the floor, not on the wall

Changes, in the sense of better function, performance, service life and technical improvements are subject to alterations.



Upgrades for Ceramill Mind

- 179151 ① Ceramill Artex (virtual articulator, contained in the pack contents of the 300 version sets)
- 179152 ② Ceramill M-Plant (Abutmentmodul)



Technical Data

Dimensions: 535 x 435 x 655mm
 Furnace chamber volume: 1l
 Electrical supply rating: 230V 50/60Hz
 Max. power consumption: 3.2kW

Ceramill Therm

178350 Ceramill Therm

Delivery volume: Ceramill Therm sintering furnace; Sintering shell, 100 x 80 x 15 mm; Sintering pearls, 200g, dia. 1mm

Accessories:

178360 ① Stackable sintering bowl, 1 pcs

Stackable only in Ceramill Therm (178350)

178361 ② Gripper for Sintering shell, stackable

178311 ③ Sintering pearls 200g, Ø 1mm



Technical Data

Dimensions: 406 x 280 x 423
 Weight: 16kg
 Electrical connection values (V/A/Hz):
 230/3,6/50-60, 115/7,0/50
 Performance: 1000W
 Suction Power: 56,6 L/sec.
 Volume: max. 52 dbA
 Filter bag size: 10 Liter
 Hose diameter: 38mm
 Hose length: 1,8 m
 Adapter diameter: conical, 37-38mm
 HEPA micro filter (97,97%),
 Filter class H12, dust class M

Ceramill Airstream

_ Suction for Ceramill Motion, Multi-x and Smartbox Invest
 _ May also be used for other devices with attached adapter

178600 Ceramill Airstream 230 V
 (100-120V = 178600V100-V120)

178610 Airstream Suction Bag (5 pcs.)

178611 Airstream Microfilter (1 pcs.)



Ceramill Scanmarker

Ceramill Scan marker is a powder spray that is applied to the surface of models or teeth to improve their visual characteristics when using a camera or scanner in the CAD/CAM technique. Suitable for extra-oral use.

760562 Ceramill Scanmarker 50ml



Ceramill Motion Starterkit

760012 Ceramill Motion Starterkit
 delivery specification: see label ★



Ceramill Motion Milling cutter

760604 Ceramill Roto Motion 1,0
 burr with magnet holder Ø 1,0mm 1 pcs. ★ 2x

760605 Ceramill Roto Motion 2,5
 burr with magnet holder Ø 2,5mm 1 pcs. ★ 2x





Ceramill Motion Blanks

760172	Ceramill ZI 71 XS, zirconium-oxide blank, dental arch shape, h=12mm	1 pcs. ★ 1x
760173	Ceramill ZI 71 S, zirconium-oxide blank, dental arch shape, h=14mm	1 pcs.
760174	① Ceramill ZI 71, zirconium-oxide blank, dental arch shape, h=16mm	1 pcs. ★ 1x
760176	Ceramill ZI 71 M, zirconium-oxide blank, dental arch shape, h=18mm	1 pcs.
760184	Ceramill ZI 71 L, zirconium-oxide blank, dental arch shape, h=20mm	1 pcs.
760175	Ceramill ZI 71 XL, zirconium-oxide blank, dental arch shape, h=25mm	1 pcs.
760301	Ceramill TEST 71 L, test blank, dental arch shape, h=20mm	1 pcs.
760307	Ceramill WAX 71 XS, wax blank, dental arch shape, h=13mm	1 pcs.
760302	② Ceramill WAX 71 L, wax blank, dental arch shape, h=20mm	1 pcs. ★ 1x
760311	Ceramill PMMA 71 XS acrylic which burns out without residue, dental arch shape, h=13mm	1 pcs.
760303	③ Ceramill PMMA 71 L, acrylic which burns out without residue, dental arch shape, h=20mm	1 pcs. ★ 1x
760309	Ceramill TEMP middle 71 XS, dyed PMMA, dental arch shape, h=13mm	1 pcs.
760305	⑤ Ceramill TEMP middle 71 L, dyed PMMA, dental arch shape, h=20mm	1 pcs. ★ 1x
760310	Ceramill TEMP light 71 XS, dyed PMMA, dental arch shape, h=13mm	1 pcs.
760306	⑥ Ceramill TEMP light 71 L, dyed PMMA, dental arch shape, h=20mm	1 pcs.
760308	Ceramill TEMP dark 71 XS, dyed PMMA dental arch shape, h=13mm	1 pcs.
760304	④ Ceramill TEMP dark 71 L, dyed PMMA, dental arch shape, h=20mm	1 pcs.



Dry solutions, size 100ml:

760471	Ceramill Liquid CL1
760472	Ceramill Liquid CL2
760473	Ceramill Liquid CL3
760474	Ceramill Liquid CL4

Ceramill Liquid

Four dilutable waterbased dye solutions for safe handling and individual colouring of the blanks.

760470	Ceramill Liquid - complete set 4 colours à 100ml, 4 colouring jar, 1 forceps	★
--------	---	---



Ceramill Marker - Excellent adhesion during processing with a water-cooled laboratory turbine. Marks imperfections without depositing anything

Ceramill Marker

Blue, oil-based contact paste for fitting zirconium oxide crowns (e.g. Ceramill ZI).

760021	Ceramill Marker, blue 3g, contact paste, 1 pcs.	★ 1x
583150	Pastebrush, brush for Ceramill Marker, 1 pcs.	★ 1x

KNOW-HOW

CERAMILL HELPDESK | TRAINING COURSES | AG WEBINARS | D-LAB 24




CERAMILL HELPDESK

The Ceramill Helpdesk: Targeted, effective help for CAD/CAM users

The Ceramill Helpdesk provides information and support from the putting into operation to the case related problem solving. A highly qualified team of dental technicians is ready to provide specific information and expert knowledge. Via phone or by desktop sharing via the Internet we are present in your laboratory and thus may offer an immediate solution that renders you productive within no time.



- _ Easy contact by ticket system - return call directly from specialists
- _ Practical support during the putting into operation
- _ Case-related problem solving
- _ Answering of questions via e-mail and phone
- _ Online support by desktop sharing - we demonstrate the next working steps on your screen
- _ Comprehensive information via Internet (www.ceramill-m-center.com)



Video tutorials on
www.ceramill-m-center.com



Highly qualified dental engineers at
Ceramill Helpdesk



Live support on your PC

TRAINING COURSES

Know the technique - master the manual skills

State-of-the-art technologies, growing expectations of patients and not least global competition step up the pace in our business sector as well. For those who want to meet the demands of the market in the long run, it is vital to be flexible and stay at the cutting edge of technology.

The entry into the digital age has also been accompanied by wide-ranging changes to training requirements and we have addressed these changes.

“Webinar” is the keyword for an efficient learning method for CAD-CAM courses that not only save you an enormous amount of effort but are mainly held on the medium which is used for working procedures. A mixture of demonstrations, training videos, live streaming and individual consultation enable you to use the system immediately after 2-3 hours.

Our contact persons are available if you have any questions or would like to register. They ensure a smooth procedure and that you can fully concentrate on your training at our locations in Germany and Austria or in the “AG training center Middle East” in Beirut (Lebanon).



Sabrina Amann



Université Antonine

INFO

AG Training Center Germany / Austria

Sabrina Amann
Intern. course organization
Fon: +43 5523 62333-57
Fax: +43 5523 55990
sabrina.amann@amangirrbach.com

For information on courses held by our dealers worldwide please contact them directly.

INFO

AG Training Center Middle East, Beirut (Lebanon)

Université Antonine
Institute of Dental Laboratory
BP 40016 - Hadath, Baabda |
Lebanon
Fon +961 1877079
Mobile +961 3133911
abdo.saleem@amangirrbach.com

TRAINING COURSES

Overview training courses Framework Management Digital

COURSE:	CAD-CAM BASIC	FMD MO1
CONTENTS:	_ System installation/set-up _ Designing on the PC (standard indications) _ Digital processing chain _ Inhouse-Outsource fabrication (incl. Motion)	
INSTRUCTOR:	AG instructors, Doris Kuster-Wohlgenannt Dental Technician	
COSTS:	EUR 390,- (1,5 days)	★★

COURSE:	CAD-CAM BASIC (Addendum Motion)	FMD MO2
CONTENTS:	_ CAM software & nesting _ Milling strategies _ Sintering process, staining, finishing	
INSTRUCTOR:	AG instructors, Doris Kuster-Wohlgenannt Dental Technician	
COSTS:	EUR 650,- (2,5 days)	★★★

COURSE:	CAD-CAM ADVANCED I	FMD MO3
CONTENTS:	_ Fabrication implant abutments _ Range of titanium bases _ Large restorations in-house digital _ Sintering process & staining	
INSTRUCTOR:	AG instructors	
COSTS:	EUR 490,- (2 days)	★★★★

COURSE:	CAD-CAM ADVANCED I	FMD MO4
CONTENTS:	_ Webinar individual course _ Ceramill M-plant abutment software basics _ Demonstration of the options for occlusally screw-retained abutments with bridges _ Entering a separate gingiva scan and scanning _ Designing of a demonstration case by the participant _ Transferring the design to the CAM Software _ Bonding the abutment with the titanium base	
INSTRUCTOR:	AG-instructors	
COSTS:	EUR 180,- (2 hours)	★★★★

COURSE:	CAD-CAM ADVANCED II	FMD MO5
CONTENTS:	_ Webinar individual course _ "individual" abutment course _ Free choice of indications by participants beforehand	
INSTRUCTOR:	AG instructors	
COSTS:	EUR 180,- (2 hours)	★★★★

COURSE:	CAD-CAM ADVANCED III	FMD MO6
CONTENTS:	_ Interface between analogue and digital dental technology in the daily routine _ Complete procedure from the model to finished coping _ Focus: The virtual articulator and all its possibilities _ Extended analysis and measuring options	
INSTRUCTOR:	AG instructors	
COSTS:		★★★★

COURSE:	CAD-CAM ADVANCED IV ZA/ZT	FMD MO7
CONTENTS:	_ Interdisciplinary dentist-dental technician _ From framework to the overall aesthetic and functional outcome _ Zirconia based	
INSTRUCTOR:	AG instructors	
COSTS:		★★★★



COURSE:	CAD-CAM ADVANCED V	FMD MO8
CONTENTS:	_ Telescope crown technique _ Fabrication of an acrylic primary crown _ Fitting primary pattern (telescope set) _ Fabrication of the secondary crown _ Fitting, feedback...	
INSTRUCTOR:	AG instructors	
COSTS:		★★★★★

Scope of services

As a matter of principle, catering for whole-day-courses or several-day-courses is included in the participation fee. The courses are acknowledged with a certificate.

Cancellation guarantee

Cancellation shall be made in written form and at least 4 weeks prior to the beginning of the event. After this period, 50% of the participation fee is due; from 1 week prior to the beginning of the course 100% of the participation fee is due. Optionally, the participation fee can be transferred to another course in the following 12 months. The participation can be assigned to a substitute free of charge. If AmannGirrbach is forced to cancel the event due to organisational or other reasons, the already paid participation fee will be reimbursed. Claims beyond this are excluded.

- ★ = Degree of difficulty (1★ = easy)
-  = Webinar
-  = Course in planning, current information: Course Hotline





powered by




netviewer

The new generation of online support via Internet

Webinar systems are similar in constitution and function to so-called "instant messenger programs" such as Skype, ICQ or the AOL instant messenger. Enhanced features such as the possibility to broadcast presentations or simulate a digital blackboard, make this an interactive online support service of the highest level.

Many problems that occur in day-to-day dental technical life can still only be solved from person to person. This is where Webinars offer a highly efficient alternative. Quick and uncomplicated help in the form of a support specialist, via Internet, directly into your laboratory.

There are many different types of problems that require personal assistance, where we can offer this type of support. Your only requirement is an active access.



- _ Help regarding our many systems, applications and products.
- _ No travel costs or waiting times for support technicians.
- _ Interactive online training possible for up to 30 participants.
- _ Cost-effective possibilities for training personnel.
- _ Quick help with installation of machinery.
- _ Quick help with maintenance and servicing.



INFOLINE

+43 5523 62333-0

or register directly at:
[www.amangirrbach.com/de/
 support/kursuebersicht/
 online-kurse-webinars/](http://www.amangirrbach.com/de/support/kursuebersicht/)



The dynamic E-learning Portal from AmannGirrbach - the digital "Product Specialist" for application techniques, daily updates, 24h online access.

The D-Lab24 is a platform that is now available for all users of our products and systems with the aim of transferring the know-how on the use of our systems and products simply, clearly and as directly as possible into the laboratory.

The requirements for the transfer of knowledge are diverse and therefore need different transport vehicles - for example techniques-instructions will be converted in the form of elaborate e-learning modules (step-by-step). Video sequences will also be used at challenging sections which will create optimal transparency. In addition, there are numerous know-how documents available in the D-Lab24 for downloading. Instructions for use, course information, case studies and basic dental laboratory know-how: in future you will find all the technical information in the D-Lab24 you require to use our products and systems in the practice correctly and efficiently.



- _Complex e-learning modules on AG methods describe every single move.
- _Processing assurance for your laboratory with the option of breaking down each single working step into the smallest details.
- _Easy, cost-effective and efficient training for your laboratory staff.
- _Videos enable complex working steps to be understood easily.
- _Downloads available on working, application and operating instructions guarantee the correct application of our products and systems from the start.
- _The entire know-how on products, systems and applications -accessible 24h
- _Always up-to-date due to continual updates
- _Extensive know-how pool - continually enhanced with new contents.



ACCESS CODE

Available forthwith from your specialist dealer, via hotline **+49 7231 957-100** Or on order via email: **d-lab24@amanngirrbach**

No online registration necessary.



SERVICE

CUSTOMER SERVICE | TECHNICAL SERVICE | SERVICE PARTNER

AG LIVE LABS

GENERAL INFORMATION



AmannGirrbach Products can only be acquired outside of Germany and Austria from our authorised dealers

You can find a specialist dealer in your area through our worldwide network of AmannGirrbach dealers.

If you do not have an active Internet access or if you cannot find a dealer in your area, please contact our export department at the AmannGirrbach headquarters in Koblach, Austria:

INFOLINE

Dealer directory:


www.amanngirrbach.com/en/contact/find-dealers/

AmannGirrbach Headquarter:

AmannGirrbach AG

Fon: +43 5523 62333-0 | Fax: +43 5523 55990

E-mail: austria@amanngirrbach.com



Authorised AmannGirrbach service partners in your area

The authorised service companies from AmannGirrbach provide a repair service using original AmannGirrbach machines and spare parts. These service companies have trained personnel and the technical information necessary for the correct and professional repair of our machines.

If you cannot find your service company, please contact our in-house technical service department:

INFOLINE

Technical Service/Servicepartner:

Hr. Jürgen Nachbaur

Fon: +43 5523 62333-207

juergen.nachbaur@amanngirrbach.com

or

www.amanngirrbach.com/en/contact/technical-service-int/



AG LIVE LABS

Experience the practical use of equipment, materials and systems, at a location near you



Information days are being held in certified laboratories worldwide. They are among the leading dental laboratories in their region in terms of size, quality, equipment and innovative capability.



Exchange of experiences between colleagues - this describes the service provided by the new AmannGirrbach Live Labs.

Registrations and current dates are available from your dealer or you can find a current list of all AG Live Labs at www.amanngirrbach.com

GENERAL INFORMATION

In this catalogue, the product information is restricted to some elementary features and applications.

You will find more information in the according system brochures and/or leaflets.

**Technical data**

In general, all dimensions are indicated in millimeters (mm) lined up in depth/length x width x height:

1. depth
2. width
3. height

If this is not the case, the differing indications are clearly defined (e.g. diameter = Ø, drilling, thickness, etc.)

Delivery volume

Introductory kits of systems or basic equipment for units are presented, on one hand, as complete set and, on the other hand, as single articles with different article numbers for repeat orders.

Replacement parts

... will be mentioned only in cases where there is high customer demand. All replacement parts are listed with article number in the unit's manual.

Changes, in the sense of better function, performance, service life and technical improvements are subject to alterations.



AMANNGIRRBACH

Amann Girrbach AG

Herrschaftswiesen 1
6842 Koblach | Austria

Fon +43 5523 62333-399
Fax +43 5523 55990

austria@amanngirrbach.com
www.amanngirrbach.com



AMANNGIRRBACH

Amann Girschbach AG

Herrschaftswiesen 1
6842 Koblach | Austria

Fon +43 5523 62333-399
Fax +43 5523 55990

austria@amanngirschbach.com
www.amanngirschbach.com



AMANNGIRRBACH

Amann Girschbach AG

Herrschaftswiesen 1
6842 Koblach | Austria

Fon +43 5523 62333-399
Fax +43 5523 55990

austria@amanngirschbach.com
www.amanngirschbach.com



AMANNGIRRBACH

Headquarter:
Amann Girrbach AG
Herrschaftswiesen 1
6842 Koblach | Austria
Fon +43 5523 62333-0
Fax +43 5523 55990

Amann Girrbach America, Inc.
12169 Villa Road
Spring Hill | FL 34609
Fon 1-800 851 3719
Fax 352 398 1409

austria@amanngirrbach.com
america@amanngirrbach.com
www.amanngirrbach.com